

cellTRAY[®] Dish by Nanopoint

The cellTRAY[®] Dish, has been designed to contain Nanopoint's cellTRAY, and any standard 1" x 3" microscope slide, during the cell loading, incubation, washing, and staining stages of an experiment.



Optically clear

Autoclavable & reusable

Frosted rim for documentation

Pedestals enable easy removal of slides

Volume similar to Petri dish

The two pedestal holders in the bottom keep the cellTRAY secure while the dish is tipped for the total extraction of media, PBS and other treatments. The lid contains a frosted area that can be used to write notes on the contents of the dish. The cellTRAY Dish is autoclavable to ensure sterility prior to use and can be stored safely in a commercial freezer without any signs of warping.

Nanopoint, Inc.
900 Fort Street Mall, Suite A20
Honolulu, HI 96813
Phone: (808) 457-1145

Order online at:
www.nanopointimaging.com

